

APPENDIX

Tabell S1. Allelvarianter inkludert i Ery Q® kit (BAG Diagnostics, Lich, Germany) for blodtypesystemene Rh, Kell, Kidd, Duffy, og MNS. Reprodusert fra intern, upublisert dokumentasjon fra BAG Diagnostics (2022), med tillatelse.

Ery Q® kit	Allel
RH	RHD*01
	RHD*01N.01 (dd)
	RHD*DEL1 (K409K)
	RHD*11 (M295I)
	RHD*DEL8 (IVS3+1G>A)
	RHD*08N.01 (Psi/ψ)
	RHD*01N.08 (W16X)
	RHD-CE (8-9)-D
	RHCE*C
	RHCE*C ^w
	RHCE*E
	RHCE*e
	RHCE*c
Partiell D	RHD*01
	RHD*01N.01
	RHD*02
	RHD*03.01,03.03
	RHD*03.02
	RHD*04.01
	RHD*04.06, *04.05
	RHD*04.03
	RHD*04.04
	RHD*05.01-*05.10
	RHD*13.01, *13.02
	RHD*06.01
	RHD*06.02
	RHD*06.03
	RHD*06.04
	RHD*07.01, 07.02
	RHD*09.01.00
	RHD*09.02.00, *09.02.01
	RHD*10.00-*10.03, *10.06 -
	*10.15
	RHD*10.04, *10.05., *10.05.01
	RHD*14.01
	RHD*14.02
	RHD*17.01 - *17.03
	RHD*17.05
	RHD*08N.01
	RHD*19
	RHD*25
	RHCE*01.22
	RHD*01N.02
	RHD*01N.03
	RHD*01N.04
	RHD*01N.05
	RHD*01N.07
	RHD*D-CE(8-9)-D
	RHD(delEx9)
	RHD*D-CE(10)
	RHD*01N.09
Svak D	RHD*01W.1.1
	RHD*01W.1
	RHD*01W.2
	RHD*01W.3
	RHD*01W.38
	RHD*01W.5
	RHD*01W.17
	RHD*01W.20
	RHD*11
	RHD*15
	RHD*01W.31
	RHD*09.01.00

*RHD*01EL.01*
*RHD*01EL.08*
*RHD*09.05*
*RHD*08N.01*
*RHD*01W.14*
*RHD*09.03.01,09.04*

KKD/MNS	<i>KEL*01.01</i>
	<i>KEL*02</i>
	<i>JK*01</i>
	<i>JK*02</i>
	<i>JK*01N.06</i>
	<i>JK*02N.06</i>
	<i>FY*01</i>
	<i>FY*02</i>
	<i>FY*02W.01 / FY*01W.01)</i> ¹
	<i>FY*02N.01 / FY*01N.01)</i> ¹
	<i>GYP A*01</i>
	<i>GYP A*02</i>
	<i>GYP B*03</i>
	<i>GYP B*04</i>
	<i>GYP B*03N.01</i>
	<i>GYP B*03N.04</i>

¹ Ved enkelte kombinasjoner av nærliggende SNP'er kan Ery Q® ikke skille mellom alternative haplotyper. Mer enn én underliggende allelkombinasjon kan derfor forklare den rapporterte genotypen.